

Antibodies Product Datasheet

Product Name:	JMY Antibody
Catalog Number:	C10003
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JMY. AA range:931-980
Concentration:	1 mg/ml
Gene Name:	JMY
Protein Name:	Junction-mediating and -regulatory protein
Other Name:	JMY; Junction-mediating and -regulatory protein
Human Swiss Prot No.:	Q8N9B5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JMY Polyclonal Antibody detects endogenous levels of JMY protein.
Storage and Stability:	-20°C/1 year