

Antibodies Product Datasheet

Product Name:	DR4 Antibody
Catalog Number:	C10014
Reactivity:	Human, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNFRSF10A. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	TNFRSF10A
Protein Name:	Tumor necrosis factor receptor superfamily member 10A
Other Name:	TNFRSF10A; APO2; DR4; TRAILR1; Tumor necrosis factor receptor superfamily member 10A; Death receptor 4; TNF-related apoptosis-inducing ligand receptor 1; TRAIL receptor 1; TRAIL-R1; CD antigen CD261
Human Swiss Prot No.:	O00220
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DR4 Polyclonal Antibody detects endogenous levels of DR4 protein.
Storage and Stability:	-20°C/1 year