

Antibodies Product Datasheet

Product Name:	FoxE1 Antibody
Catalog Number:	C10022
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TTF2. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	FOXE1
Protein Name:	Forkhead box protein E1
Other Name:	FOXE1; FKHL15; FOXE2; TTF2; TTF2; Forkhead box protein E1; Forkhead box protein E2; Forkhead-related protein FKHL15; HFKH4; HNF-3/fork head-like protein 5; HFKL5; Thyroid transcription factor 2; TTF-2
Human Swiss Prot No.:	O00358
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxE1 Polyclonal Antibody detects endogenous levels of FoxE1 protein.
Storage and Stability:	-20°C/1 year