

Antibodies Product Datasheet

Product Name:	DOC-1 Antibody
Catalog Number:	C10037
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDKAP1. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	CDK2AP1
Protein Name:	Cyclin-dependent kinase 2-associated protein 1
Other Name:	CDK2AP1; CDKAP1; DOC1; Cyclin-dependent kinase 2-associated protein 1; CDK2-associated protein 1; Deleted in oral cancer 1; DOC-1; Putative oral cancer suppressor
Human Swiss Prot No.:	O14519
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DOC-1 Polyclonal Antibody detects endogenous levels of DOC-1 protein.
Storage and Stability:	-20°C/1 year