

Antibodies Product Datasheet

Product Name:	Hec1 Antibody
Catalog Number:	C10051
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KNTC2. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	NDC80
Protein Name:	Kinetochore protein NDC80 homolog
Other Name:	NDC80; HEC; HEC1; KNTC2; Kinetochore protein NDC80 homolog; Highly expressed in cancer protein; Kinetochore protein Hec1; HsHec1; Kinetochore-associated protein 2; Retinoblastoma-associated protein HEC
Human Swiss Prot No.:	O14777
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hec1 Polyclonal Antibody detects endogenous levels of Hec1 protein.
Storage and Stability:	-20°C/1 year