

Antibodies Product Datasheet

Product Name:	LIGHT Antibody
Catalog Number:	C10109
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNF14. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	TNFSF14
Protein Name:	Tumor necrosis factor ligand superfamily member 14
Other Name:	TNFSF14; HVEM-L; LIGHT; Tumor necrosis factor ligand superfamily member 14; Herpes virus entry mediator ligand; HVEM-L; Herpesvirus entry mediator ligand; CD antigen CD258
Human Swiss Prot No.:	O43557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LIGHT Polyclonal Antibody detects endogenous levels of LIGHT protein.
Storage and Stability:	-20°C/1 year