

Antibodies Product Datasheet

Product Name:	CRSP150 Antibody
Catalog Number:	C10125
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED14. AA range:701-750
Concentration:	1 mg/ml
Gene Name:	MED14
Protein Name:	Mediator of RNA polymerase II transcription subunit 14
Other Name:	MED14; ARC150; CRSP2; CXorf4; DRIP150; EXLM1; RGR1; TRAP170; Mediator of RNA polymerase II transcription subunit 14; Activator-recruited cofactor 150 kDa component; ARC150; Cofactor required for Sp1 transcriptional activation subunit 2; CRS
Human Swiss Prot No.:	O60244
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRSP150 Polyclonal Antibody detects endogenous levels of CRSP150 protein.
Storage and Stability:	-20°C/1 year