

Antibodies Product Datasheet

Product Name:	BubR1 Antibody
Catalog Number:	C10142
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BUB1B. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	BUB1B
Protein Name:	Mitotic checkpoint serine/threonine-protein kinase BUB1 beta
Other Name:	BUB1B; BUBR1; MAD3L; SSK1; Mitotic checkpoint serine/threonine-protein kinase BUB1 beta; MAD3/BUB1-related protein kinase; hBUBR1; Mitotic checkpoint kinase MAD3L; Protein SSK1
Human Swiss Prot No.:	O60566
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BubR1 Polyclonal Antibody detects endogenous levels of BubR1 protein.
Storage and Stability:	-20°C/1 year