

Antibodies Product Datasheet

Product Name:	DNA pol ζ Antibody
Catalog Number:	C10144
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DNA Polymerase zeta. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	REV3L
Protein Name:	DNA polymerase zeta catalytic subunit
Other Name:	REV3L; POLZ; REV3; DNA polymerase zeta catalytic subunit; Protein reversionless 3-like; REV3-like; hREV3
Human Swiss Prot No.:	O60673
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DNA pol ζ Polyclonal Antibody detects endogenous levels of DNA pol ζ protein.
Storage and Stability:	-20°C/1 year