

Antibodies Product Datasheet

Product Name:	SLC43A1 Antibody
Catalog Number:	C10170
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAT3. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	SLC43A1
Protein Name:	Large neutral amino acids transporter small subunit 3
Other Name:	SLC43A1; LAT3; PB39; POV1; Large neutral amino acids transporter small subunit 3; L-type amino acid transporter 3; Prostate cancer overexpressed gene 1 protein; Solute carrier family 43 member 1
Human Swiss Prot No.:	O75387
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLC43A1 Polyclonal Antibody detects endogenous levels of SLC43A1 protein.
Storage and Stability:	-20°C/1 year