

Antibodies Product Datasheet

Product Name:	FAST-1/2 Antibody
Catalog Number:	C10188
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXH1. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	FOXH1
Protein Name:	Forkhead box protein H1
Other Name:	FOXH1; FAST1; FAST2; Forkhead box protein H1; Forkhead activin signal transducer 1; Fast-1; hFAST-1; Forkhead activin signal transducer 2; Fast-2
Human Swiss Prot No.:	O75593
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAST-1/2 Polyclonal Antibody detects endogenous levels of FAST-1/2 protein.
Storage and Stability:	-20°C/1 year