

Antibodies Product Datasheet

Product Name:	DOC-1R Antibody
Catalog Number:	C10199
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDKA2. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	CDK2AP2
Protein Name:	Cyclin-dependent kinase 2-associated protein 2
Other Name:	CDK2AP2; DOC1R; Cyclin-dependent kinase 2-associated protein 2; CDK2-associated protein 2; DOC-1-related protein; DOC-1R
Human Swiss Prot No.:	O75956
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DOC-1R Polyclonal Antibody detects endogenous levels of DOC-1R protein.
Storage and Stability:	-20°C/1 year