

Antibodies Product Datasheet

Product Name:	AKAP 110 Antibody
Catalog Number:	C10200
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKAP3. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	AKAP3
Protein Name:	A-kinase anchor protein 3
Other Name:	AKAP3; AKAP110; SOB1; A-kinase anchor protein 3; AKAP-3; A-kinase anchor protein 110 kDa; AKAP 110; Cancer/testis antigen 82; CT82; Fibrous sheath protein of 95 kDa; FSP95; Fibrousheathin I; Fibrousheathin-1; Protein kinase A-anchoring prot
Human Swiss Prot No.:	O75969
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AKAP 110 Polyclonal Antibody detects endogenous levels of AKAP 110 protein.
Storage and Stability:	-20°C/1 year