

Antibodies Product Datasheet

Product Name:	MEK Kinase-6 Antibody
Catalog Number:	C10227
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAP3K6. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	MAP3K6
Protein Name:	Mitogen-activated protein kinase kinase kinase 6
Other Name:	MAP3K6; ASK2; MAPKKK6; MEKK6; Mitogen-activated protein kinase kinase kinase 6; Apoptosis signal-regulating kinase 2
Human Swiss Prot No.:	O95382
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEK Kinase-6 Polyclonal Antibody detects endogenous levels of MEK Kinase-6 protein.
Storage and Stability:	-20°C/1 year