

Antibodies Product Datasheet

Product Name:	KLF8 Antibody
Catalog Number:	C10236
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KLF8. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	KLF8
Protein Name:	Krueppel-like factor 8
Other Name:	KLF8; BKLF3; ZNF741; Krueppel-like factor 8; Basic krueppel-like factor 3; Zinc finger protein 741
Human Swiss Prot No.:	O95600
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KLF8 Polyclonal Antibody detects endogenous levels of KLF8 protein.
Storage and Stability:	-20°C/1 year