

Antibodies Product Datasheet

Product Name:	HGK Antibody
Catalog Number:	C10242
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEKKK 4. AA range:406-455
Concentration:	1 mg/ml
Gene Name:	MAP4K4
Protein Name:	Mitogen-activated protein kinase kinase kinase 4
Other Name:	MAP4K4; HGK; KIAA0687; NIK; Mitogen-activated protein kinase kinase kinase 4; HPK/GCK-like kinase HGK; MAPK/ERK kinase kinase 4; MEK kinase kinase 4; MEKKK 4; Nck-interacting kinase
Human Swiss Prot No.:	O95819
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HGK Polyclonal Antibody detects endogenous levels of HGK protein.
Storage and Stability:	-20°C/1 year