

Antibodies Product Datasheet

Product Name:	HAND1 Antibody
Catalog Number:	C10253
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HAND1. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	HAND1
Protein Name:	Heart- and neural crest derivatives-expressed protein 1
Other Name:	HAND1; BHLHA27; EHAND; Heart- and neural crest derivatives-expressed protein 1; Class A basic helix-loop-helix protein 27; bHLHa27; Extraembryonic tissues; heart, autonomic nervous system and neural crest derivatives-expressed protein 1; eH
Human Swiss Prot No.:	O96004
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HAND1 Polyclonal Antibody detects endogenous levels of HAND1 protein.
Storage and Stability:	-20°C/1 year