

Antibodies Product Datasheet

Product Name:	HC-II Antibody
Catalog Number:	C10281
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Heparin Cofactor II. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	SERPIND1
Protein Name:	Heparin cofactor 2
Other Name:	SERPIND1; HCF2; Heparin cofactor 2; Heparin cofactor II; HC-II; Protease inhibitor leuserpin-2; HLS2; Serpin D1
Human Swiss Prot No.:	P05546
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HC-II Polyclonal Antibody detects endogenous levels of HC-II protein.
Storage and Stability:	-20°C/1 year