

Antibodies Product Datasheet

Product Name:	Hox-C6 Antibody
Catalog Number:	C10318
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HOXC6. AA range:186-235
Concentration:	1 mg/ml
Gene Name:	HOXC6
Protein Name:	Homeobox protein Hox-C6
Other Name:	HOXC6; HOX3C; Homeobox protein Hox-C6; Homeobox protein CP25; Homeobox protein HHO.C8; Homeobox protein Hox-3C
Human Swiss Prot No.:	P09630
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hox-C6 Polyclonal Antibody detects endogenous levels of Hox-C6 protein.
Storage and Stability:	-20°C/1 year