

Antibodies Product Datasheet

Product Name:	HKR1 Antibody
Catalog Number:	C10329
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HKR1. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	HKR1
Protein Name:	Krueppel-related zinc finger protein 1
Other Name:	HKR1; ZNF875; Krueppel-related zinc finger protein 1; Protein HKR1; Zinc finger protein 875
Human Swiss Prot No.:	P10072
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HKR1 Polyclonal Antibody detects endogenous levels of HKR1 protein.
Storage and Stability:	-20°C/1 year