

Antibodies Product Datasheet

Product Name:	A-Myb Antibody
Catalog Number:	C10332
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYB-A. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	MYBL1
Protein Name:	Myb-related protein A
Other Name:	MYBL1; AMYB; Myb-related protein A; A-Myb; Myb-like protein 1
Human Swiss Prot No.:	P10243
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A-Myb Polyclonal Antibody detects endogenous levels of A-Myb protein.
Storage and Stability:	-20°C/1 year