

Antibodies Product Datasheet

Product Name:	KOX11 Antibody
Catalog Number:	C10394
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF18. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	ZNF18
Protein Name:	Zinc finger protein 18
Other Name:	ZNF18; HDSG1; KOX11; ZKSCAN6; ZNF535; Zinc finger protein 18; Heart development-specific gene 1 protein; Zinc finger protein 535; Zinc finger protein KOX11; Zinc finger protein with KRAB and SCAN domains 6
Human Swiss Prot No.:	P17022
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KOX11 Polyclonal Antibody detects endogenous levels of KOX11 protein.
Storage and Stability:	-20°C/1 year