

Antibodies Product Datasheet

Product Name:	KOX17 Antibody
Catalog Number:	C10397
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF24. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	ZNF24
Protein Name:	Zinc finger protein 24
Other Name:	ZNF24; KOX17; ZNF191; ZSCAN3; Zinc finger protein 24; Retinoic acid suppression protein A; RSG-A; Zinc finger and SCAN domain-containing protein 3; Zinc finger protein 191; Zinc finger protein KOX17
Human Swiss Prot No.:	P17028
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KOX17 Polyclonal Antibody detects endogenous levels of KOX17 protein.
Storage and Stability:	-20°C/1 year