

Antibodies Product Datasheet

Product Name:	Egr-2 Antibody
Catalog Number:	C10416
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EGR1/2. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	EGR2
Protein Name:	Early growth response protein 2
Other Name:	EGR1; KROX24; ZNF225; Early growth response protein 1; EGR-1; AT225; Nerve growth factor-induced protein A; NGFI-A; Transcription factor ETR103; Transcription factor Zif268; Zinc finger protein 225; Zinc finger protein Krox-24; EGR2; KROX20
Human Swiss Prot No.:	P18146/P11161
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Egr-2 Polyclonal Antibody detects endogenous levels of Egr-2 protein.
Storage and Stability:	-20°C/1 year