

Antibodies Product Datasheet

Product Name:	Ephrin-A1 Antibody
Catalog Number:	C10440
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EFNA1. AA range:66-115
Concentration:	1 mg/ml
Gene Name:	EFNA1
Protein Name:	Ephrin-A1
Other Name:	EFNA1; EPLG1; LERK1; TNFAIP4; Ephrin-A1; EPH-related receptor tyrosine kinase ligand 1; LERK-1; Immediate early response protein B61; Tumor necrosis factor alpha-induced protein 4; TNF alpha-induced protein 4
Human Swiss Prot No.:	P20827
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ephrin-A1 Polyclonal Antibody detects endogenous levels of Ephrin-A1 protein.
Storage and Stability:	-20°C/1 year