

Antibodies Product Datasheet

Product Name:	A20 Antibody
Catalog Number:	C10444
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNAP3. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	TNFAIP3
Protein Name:	Tumor necrosis factor alpha-induced protein 3
Other Name:	TNFAIP3; OTUD7C; Tumor necrosis factor alpha-induced protein 3; TNF alpha-induced protein 3; OTU domain-containing protein 7C; Putative DNA-binding protein A20; Zinc finger protein A20
Human Swiss Prot No.:	P21580
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A20 Polyclonal Antibody detects endogenous levels of A20 protein.
Storage and Stability:	-20°C/1 year