

Antibodies Product Datasheet

Product Name:	DGK- α Antibody
Catalog Number:	C10463
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKA. AA range:304-353
Concentration:	1 mg/ml
Gene Name:	DGKA
Protein Name:	Diacylglycerol kinase alpha
Other Name:	DGKA; DAGK; DAGK1; Diacylglycerol kinase alpha; DAG kinase alpha; 80 kDa diacylglycerol kinase; Diglyceride kinase alpha; DGK-alpha
Human Swiss Prot No.:	P23743
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- α Polyclonal Antibody detects endogenous levels of DGK- α protein.
Storage and Stability:	-20°C/1 year