

Antibodies Product Datasheet

Product Name:	Cdk2 Antibody
Catalog Number:	C10474
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDK2. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	CDK2
Protein Name:	Cyclin-dependent kinase 2
Other Name:	CDK2; CDKN2; Cyclin-dependent kinase 2; Cell division protein kinase 2; p33 protein kinase
Human Swiss Prot No.:	P24941
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdk2 Polyclonal Antibody detects endogenous levels of Cdk2 protein.
Storage and Stability:	-20°C/1 year