

Antibodies Product Datasheet

Product Name:	HuB/D Antibody
Catalog Number:	C10486
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELAV2/4. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	ELAVL2/ELAVL4
Protein Name:	ELAV-like protein 2/ELAV-like protein 4
Other Name:	ELAVL2; HUB; ELAV-like protein 2; ELAV-like neuronal protein 1; Hu-antigen B; HuB; Nervous system-specific RNA-binding protein Hel-N1; ELAVL4; HUD; PNEM; ELAV-like protein 4; Hu-antigen D; HuD; Paraneoplastic encephalomyelitis antigen HuD
Human Swiss Prot No.:	Q12926/P26378
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HuB/D Polyclonal Antibody detects endogenous levels of HuB/D protein.
Storage and Stability:	-20°C/1 year