

Antibodies Product Datasheet

Product Name:	Cyclophilin F Antibody
Catalog Number:	C10516
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPIF. AA range:86-135
Concentration:	1 mg/ml
Gene Name:	PPIF
Protein Name:	Peptidyl-prolyl cis-trans isomerase F mitochondrial
Other Name:	PPIF; CYP3; Peptidyl-prolyl cis-trans isomerase F; mitochondrial; PPIase F; Cyclophilin F; Rotamase F
Human Swiss Prot No.:	P30405
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclophilin F Polyclonal Antibody detects endogenous levels of Cyclophilin F protein.
Storage and Stability:	-20°C/1 year