

Antibodies Product Datasheet

Product Name:	NK-TR Antibody
Catalog Number:	C10518
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NKTR. AA range:784-833
Concentration:	1 mg/ml
Gene Name:	NKTR
Protein Name:	NK-tumor recognition protein
Other Name:	NKTR; NK-tumor recognition protein; NK-TR protein; Natural-killer cells cyclophilin-related protein
Human Swiss Prot No.:	P30414
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NK-TR Polyclonal Antibody detects endogenous levels of NK-TR protein.
Storage and Stability:	-20°C/1 year