

Antibodies Product Datasheet

Product Name:	Coronin 1A Antibody
Catalog Number:	C10520
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Coronin 1A . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	CORO1A
Protein Name:	Coronin-1A
Other Name:	CORO1A; CORO1; Coronin-1A; Coronin-like protein A; Clipin-A; Coronin-like protein p57; Tryptophan aspartate-containing coat protein; TACO
Human Swiss Prot No.:	P31146
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Coronin 1A Polyclonal Antibody detects endogenous levels of Coronin 1A protein.
Storage and Stability:	-20°C/1 year