

Antibodies Product Datasheet

Product Name:	GRK 4 Antibody
Catalog Number:	C10536
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from GRK 4 . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	GRK4
Protein Name:	G protein-coupled receptor kinase 4
Other Name:	GRK4; GPRK2L; GPRK4; G protein-coupled receptor kinase 4; G protein-coupled receptor kinase GRK4; ITI1
Human Swiss Prot No.:	P32298
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRK 4 Polyclonal Antibody detects endogenous levels of GRK 4 protein.
Storage and Stability:	-20°C/1 year