

Antibodies Product Datasheet

Product Name:	CD30-L Antibody
Catalog Number:	C10540
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD153. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	TNFSF8
Protein Name:	Tumor necrosis factor ligand superfamily member 8
Other Name:	TNFSF8; CD30L; CD30LG; Tumor necrosis factor ligand superfamily member 8; CD30 ligand; CD30-L; CD antigen CD153
Human Swiss Prot No.:	P32971
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD30-L Polyclonal Antibody detects endogenous levels of CD30-L protein.
Storage and Stability:	-20°C/1 year