

Antibodies Product Datasheet

Product Name:	CstF-64 Antibody
Catalog Number:	C10542
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSTF2. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	CSTF2
Protein Name:	Cleavage stimulation factor subunit 2
Other Name:	CSTF2; Cleavage stimulation factor subunit 2; CF-1 64 kDa subunit; Cleavage stimulation factor 64 kDa subunit; CSTF 64 kDa subunit; CstF-64
Human Swiss Prot No.:	P33240
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CstF-64 Polyclonal Antibody detects endogenous levels of CstF-64 protein.
Storage and Stability:	-20°C/1 year