

Antibodies Product Datasheet

Product Name:	DAN Antibody
Catalog Number:	C10599
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NBL1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	NBL1
Protein Name:	Neuroblastoma suppressor of tumorigenicity 1
Other Name:	NBL1; DAN; DAND1; Neuroblastoma suppressor of tumorigenicity 1; DAN domain family member 1; Protein N03; Zinc finger protein DAN
Human Swiss Prot No.:	P41271
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DAN Polyclonal Antibody detects endogenous levels of DAN protein.
Storage and Stability:	-20°C/1 year