

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | MAGE-9 Antibody                                                                                                       |
| <b>Catalog Number:</b>        | C10615                                                                                                                |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                     |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA                                                                                    |
| <b>Dilutions:</b>             | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.                            |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from MAGE-9 . at AA range: 1-80                                                           |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | MAGEA9                                                                                                                |
| <b>Protein Name:</b>          | Melanoma-associated antigen 9                                                                                         |
| <b>Other Name:</b>            | MAGEA9; MAGE9; MAGEA9A; MAGEA9B; Melanoma-associated antigen 9; Cancer/testis antigen 1.9; CT1.9; MAGE-9 antigen      |
| <b>Human Swiss Prot No.:</b>  | P43362                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | MAGE-9 Polyclonal Antibody detects endogenous levels of MAGE-9 protein.                                               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |