

Antibodies Product Datasheet

Product Name:	CLK1 Antibody
Catalog Number:	C10651
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLK1. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	CLK1
Protein Name:	Dual specificity protein kinase CLK1
Other Name:	CLK1; CLK; Dual specificity protein kinase CLK1; CDC-like kinase 1
Human Swiss Prot No.:	P49759
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLK1 Polyclonal Antibody detects endogenous levels of CLK1 protein.
Storage and Stability:	-20°C/1 year