

Antibodies Product Datasheet

Product Name:	Hint1 Antibody
Catalog Number:	C10655
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HINT1. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	HINT1
Protein Name:	Histidine triad nucleotide-binding protein 1
Other Name:	HINT1; HINT; PKCI1; PRKCNH1; Histidine triad nucleotide-binding protein 1; Adenosine 5'-monophosphoramidase; Protein kinase C inhibitor 1; Protein kinase C-interacting protein 1; PKCI-1
Human Swiss Prot No.:	P49773
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hint1 Polyclonal Antibody detects endogenous levels of Hint1 protein.
Storage and Stability:	-20°C/1 year