

Antibodies Product Datasheet

Product Name:	ROR γ Antibody
Catalog Number:	C10676
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RORG. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	RORC
Protein Name:	Nuclear receptor ROR-gamma
Other Name:	RORC; NR1F3; RORG; RZRG; Nuclear receptor ROR-gamma; Nuclear receptor RZR-gamma; Nuclear receptor subfamily 1 group F member 3; Retinoid-related orphan receptor-gamma
Human Swiss Prot No.:	P51449
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ROR γ Polyclonal Antibody detects endogenous levels of ROR γ protein.
Storage and Stability:	-20°C/1 year