

Antibodies Product Datasheet

Product Name:	IRX3 Antibody
Catalog Number:	C10796
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IRX3. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	IRX3
Protein Name:	Iroquois-class homeodomain protein IRX-3
Other Name:	IRX3; IRXB1; Iroquois-class homeodomain protein IRX-3; Homeodomain protein IRXB1; Iroquois homeobox protein 3
Human Swiss Prot No.:	P78415
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IRX3 Polyclonal Antibody detects endogenous levels of IRX3 protein.
Storage and Stability:	-20°C/1 year