

Antibodies Product Datasheet

Product Name:	CRBP-III Antibody
Catalog Number:	C10805
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CRBP III. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	RBP5
Protein Name:	Retinol-binding protein 5
Other Name:	RBP5; Retinol-binding protein 5; Cellular retinol-binding protein III; CRBP-III; HRBPiso
Human Swiss Prot No.:	P82980
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRBP-III Polyclonal Antibody detects endogenous levels of CRBP-III protein.
Storage and Stability:	-20°C/1 year