

Antibodies Product Datasheet

Product Name:	mtTFA Antibody
Catalog Number:	C10817
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TFAM. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	TFAM
Protein Name:	Transcription factor A mitochondrial
Other Name:	TFAM; TCF6; TCF6L2; Transcription factor A; mitochondrial; mtTFA; Mitochondrial transcription factor 1; MtTF1; Transcription factor 6; TCF-6; Transcription factor 6-like 2
Human Swiss Prot No.:	Q00059
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	mtTFA Polyclonal Antibody detects endogenous levels of mtTFA protein.
Storage and Stability:	-20°C/1 year