

Antibodies Product Datasheet

Product Name:	CstF-50 Antibody
Catalog Number:	C10882
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSTF1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CSTF1
Protein Name:	Cleavage stimulation factor subunit 1
Other Name:	CSTF1; Cleavage stimulation factor subunit 1; CF-1 50 kDa subunit; Cleavage stimulation factor 50 kDa subunit; CSTF 50 kDa subunit; CstF-50
Human Swiss Prot No.:	Q05048
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CstF-50 Polyclonal Antibody detects endogenous levels of CstF-50 protein.
Storage and Stability:	-20°C/1 year