

Antibodies Product Datasheet

Product Name:	FoxC1/2 Antibody
Catalog Number:	C10937
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXC1/2. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	FOXC1/FOXC2
Protein Name:	Forkhead box protein C1/2
Other Name:	FOXC1; FKHL7; FREAC3; Forkhead box protein C1; Forkhead-related protein FKHL7; Forkhead-related transcription factor 3; FREAC-3; FOXC2; FKHL14; MFH1; Forkhead box protein C2; Forkhead-related protein FKHL14; Mesenchyme fork head protein 1;
Human Swiss Prot No.:	Q12948/Q99958
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxC1/2 Polyclonal Antibody detects endogenous levels of FoxC1/2 protein.
Storage and Stability:	-20°C/1 year