

Antibodies Product Datasheet

Product Name:	FoxI1 Antibody
Catalog Number:	C10938
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXI1. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	FOXI1
Protein Name:	Forkhead box protein I1
Other Name:	FOXI1; FKHL10; FREAC6; Forkhead box protein I1; Forkhead-related protein FKHL10; Forkhead-related transcription factor 6; FREAC-6; Hepatocyte nuclear factor 3 forkhead homolog 3; HFH-3; HNF-3/forkhead homolog 3
Human Swiss Prot No.:	Q12951
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxI1 Polyclonal Antibody detects endogenous levels of FoxI1 protein.
Storage and Stability:	-20°C/1 year