

Antibodies Product Datasheet

Product Name:	FoxL1 Antibody
Catalog Number:	C10939
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXL1. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	FOXL1
Protein Name:	Forkhead box protein L1
Other Name:	FOXL1; FKHL11; FREAC7; Forkhead box protein L1; Forkhead-related protein FKHL11; Forkhead-related transcription factor 7; FREAC-7
Human Swiss Prot No.:	Q12952
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxL1 Polyclonal Antibody detects endogenous levels of FoxL1 protein.
Storage and Stability:	-20°C/1 year