

Antibodies Product Datasheet

Product Name:	CAF-1 p150 Antibody
Catalog Number:	C10945
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CAF-1 p150 . at AA range: 300-380
Concentration:	1 mg/ml
Gene Name:	CHAF1A
Protein Name:	Chromatin assembly factor 1 subunit A
Other Name:	CHAF1A; CAF; CAF1P150; Chromatin assembly factor 1 subunit A; CAF-1 subunit A; Chromatin assembly factor I p150 subunit; CAF-I 150 kDa subunit; CAF-I p150; hp150
Human Swiss Prot No.:	Q13111
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CAF-1 p150 Polyclonal Antibody detects endogenous levels of CAF-1 p150 protein.
Storage and Stability:	-20°C/1 year