

Antibodies Product Datasheet

Product Name:	DGK- ζ Antibody
Catalog Number:	C10999
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKZ. AA range:601-650
Concentration:	1 mg/ml
Gene Name:	DGKZ
Protein Name:	Diacylglycerol kinase zeta
Other Name:	DGKZ; DAGK6; Diacylglycerol kinase zeta; DAG kinase zeta; Diglyceride kinase zeta; DGK-zeta
Human Swiss Prot No.:	Q13574
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- ζ Polyclonal Antibody detects endogenous levels of DGK- ζ protein.
Storage and Stability:	-20°C/1 year