

Antibodies Product Datasheet

Product Name:	GK1/3 Antibody
Catalog Number:	C11031
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GK3. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	GK/GK3P
Protein Name:	Putative glycerol kinase 3
Other Name:	GK; Glycerol kinase; GK; Glycerokinase; ATP:glycerol 3-phosphotransferase; GK3P; GKP3; GKTB; Putative glycerol kinase 3; GK 3; Glycerokinase 3; ATP:glycerol 3-phosphotransferase 3; Glycerol kinase; testis specific 1
Human Swiss Prot No.:	P32189/Q14409
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GK1/3 Polyclonal Antibody detects endogenous levels of GK1/3 protein.
Storage and Stability:	-20°C/1 year